

**RACECOURSE COVID VACCINATION SITE RESUS SESSION CHECKLIST**

<b>DAILY CHECKS</b>	<b>INITIAL</b>
<b><u>Duty Resuscitation Team/Outside Responder</u></b> 1. Team members have been identified and attended process briefing	
<b><u>Check defibs (on trolley and OR box)</u></b> 1. Check that green light is flashing – if it is not this means that the AED is not ready for use - immediately alert duty manager. 2. Check pads are attached and in date. 3. Check battery is in date.	
<b><u>Check oxygen cylinders</u></b> 1. Is cylinder full? 2. Is the oxygen flowing when switched on? 3. Adult Hudson/non-rebreathe mask attached. 4. Bag contains normal adult mask+tubing and child mask+tubing.	
<b><u>Check resus trolley</u></b> 1. All contents present and within expiry date(see below). 2. Wipe down surfaces and touchpoints of trolley with Clinell wipe.	
<b><u>Check Outside Responder box</u></b> 1. All contents present and within expiry date(see below). 2. Wipe down surfaces and touchpoints with Clinell wipe.	

<b>RESUS TROLLEY CONTENTS</b>	<b>TICK</b>
X2 Sealed PPE crashpacks within expiry date	
<b>X2 Anaphylactic shock kits (sealed and in date):</b>  X2 ampoules adrenaline 1mg/ml (1:1000) X4 23G (blue) needles X4 graduated 1ml syringes	
Salbutamol inhaler	
Volumatic spacer	
Defib (see above)	
Spare adhesive defibrillator pads.	
Towel	
X2 Razors	
X2 Tuffkut shears	
Cotton wool balls	
Oxygen Cylinder (as above)	
X2 Spare Adult Hudson Masks	
X1 Spare Child Hudson Mask	
Stethoscope	
BP machine	
Pulse Oximeter	
Anaphlaxis Algorithm 2021 Laminare	
BLS Algorithm Laminare	
Clipboard with pen and X2 copies of Resus event log	

**PLEASE TURN OVER TO COMPLETE CHECKS**

OUTSIDE RESPONDER BOX CONTENTS	TICK
Gloves and aprons (black compartments in top section)	
<b>X1 Anaphylactic shock kit (sealed and in date):</b>  X3 ampoules adrenaline 1mg/ml (1:1000) X3 23G (blue) needles X3 graduated 1ml syringes Cotton wool babs	
Defib (see above)	
X1 Tuffkut shears (in defib)	
Steristrips	
Saline	
Forceps	
Sterile scissors	
Dressing packs	
Assortment of mepores and silicone dressings	
Assortment of non-adhesive wound contact dressings	
Crepe bandages	
Tape	
Stethoscope	
Sphygmomanometer	
Pulse Oximeter	
Anaphlaxis Algorithm 2021 Laminate	
BLS Algorithm Laminate	
X1 copy Resus event log +pen	

<b>NAME:</b>	<b>SIGNATURE:</b>	<b>DATE:</b>
--------------	-------------------	--------------